

The 'QUANTUM' WET CHEMISTRY ANALYSER

WHY WE ARE DIFFERENT FROM THE REST!

'QUANTUM'	OTHER ANALYSERS
✓ L.E.D technology means a longer life, stability and zero noise interference	X Use lamps that can be unstable and have a short shelf life
✓ Truly easy to use with simple processes and software ~ anyone can use it!	X Can be less user friend with complicated software and lengthy start up procedures.
✓ Ready to use reagents in most cases with an 'open' system	X May have to prepare most reagents before use
✓ Reagent stability up to 6 months (even without onboard refrigeration)	X Don't give stability even with onboard refrigeration
✓ Independent probes for both sample & reagents resulting in increased accuracy.	X Only have one probe that is used for both sampling and reagents
✓ Disposable cuvette for each test ~ never reused	X May re-wash cuvettes or use a flow cell system that can compromise tests
✓ Removable reagent rotor for refrigeration	X Are fixed so you have to individually remove reagents to the fridge and cap them
✓ 25 to 50 reagent positions ~ (subject to sample numbers)	X Have 20 - 30 positions
✓ 25 to 50 sample positions ~ (subject to reagent numbers)	X Have 10 - 20 positions
✓ Windows Operating System	X Have either DOS or designed software that is inflexible
✓ Multiple species capability ~ (20)	X Are limited to 3 - 6 species
✓ Dedicated 'Veterinary' software	X Use human 'fields'
✓ Open methodology (no fixed times)	X Have fixed times for differing methods
✓ Test capabilities only limited to what is available in the market-place	X Are limited due to fixed parameters & closed systems
✓ Networking capabilities	X Claim to network but in reality they don't!

Web: www.quantumvetdiagnostics.com Email: sales@quantumvetdiagnostics.com

Tel: 01798 813659 Fax: 01798 815273

Unit 6 & 7 • Crays Court • Crays Lane • Goose Green • Pulborough • West Sussex • RH20 2GU • United Kingdom