

CT3000i ECG Machine

Available for hire or purchase for your clinical trial.



Overview

The CT3000i ECG machine is easy to use it features user friendly controls and settings and comes complete with an Interpretation and Measurements program.

Ideal for use in a wide variety of clinical settings.

12 lead ECG with 3 channel high quality interpretive presentation.

90mm Z fold paper print out.

Specification

- Measurement Table: RR, P, PR, QRS, QT and QTC intervals, and the P, QRS and T axis
- Power: 115/240 VAC, 50/60 HZ
- Paper Speed: 5, 25, 50 mm/sec
- Battery: Built in 12 V lead acid battery, rechargeable
- Weight: 6.4 lbs/ 2.9kgs (with battery installed)
- Dimensions: 3 x 1.8 x 8.5 inches / 77.5h x 296w x 220d mm
- Sensitivity 5/10/20 mm/mV

Operating Features

- Auto-start button automatically produces the ECG traces complete with Rhythm Strip Measurements and Interpretation
- Filters for baseline, electrical and muscle tremor interference, providing clear traces
- Lead off warning and identification
- Useful copy button for extra ECG traces
- 12 simultaneously acquired ECG traces printed on 3.6 inches(90 mm) paper, eliminates trace overlap
- Paper tray accommodates Z-fold or roll paper
- Manual button provides real time print-outs
- Mains and battery operation.



Please contact us for a quotation for the hire or purchase of this product

Product description: CT3000i ECG Machine

Product code - WD7ECG001

Version: V3 10/19



Product specifications can be subject to change without prior notice

UK Office

t: +44 (0)8456 777001 f: +44 (0)8456 777002

USA Office

t: 1 800 471 9200 / 1 718 606 0516 f: 1 718 606 0955

e: hello@woodleytrialsolutions.com

www.woodleytrialsolutions.com

woodley
TRIAL SOLUTIONS

Woodley Trial Solutions is a division of Woodley Equipment Company